

# CO2

## ALTERNATE NAMES

## SPECIMEN TYPE

Serum or heparin plasma

## STORAGE

Refrigerate

## ROUTINE / STAT T.A.T.

24 hours / 30 mins

## REFERENCE RANGE

0-28 days: 13-22 mmol/L  
29 days – 17 yrs: 20-28 mmol/L  
18 yrs – 60 yrs: 22-29 mmol/L  
>60 yrs: 23-31 mmol/L

## STABILITY

1 week refrigerated

## SPECIAL INSTRUCTIONS