

INFLUENZA/RSV PCR

ALTERNATE NAMES

SPECIMEN TYPE

NASAL ASPIRATES /WASHES OR
NASOPHARYNEAL SWABS
COLLECTED AND PLACED IN
UNIVERSAL TRANSPORT MEDIUM.

STORAGE

2-8°C. SAMPLES MAY BE STORED UP
TO 72 HOURS BEFORE TESTING

ROUTINE / STAT T.A.T.

2 HOURS

REFERENCE RANGE

INFLUENZA A: NEGATIVE
INFLUENZA B: NEGATIVE
RSV: NEGATIVE

STABILITY

SPECIAL INSTRUCTIONS