

Crystal Identification

ALTERNATE NAMES

Crystal ID, Uric Crystal ID,
Monosodium urate crystals

SPECIMEN TYPE

Synovial fluid in a sterile container
or a red top tube

STORAGE / STABILITY

8 hours Room Temp
24 hours Refrigerated

ROUTINE / STAT T.A.T.

24 hours / 60 minutes

REFERENCE RANGE

Negative

**The presence or absence of both Uric Acid (Monosodium Urate) crystals and CPPD (calcium pyrophosphate dehydrate) crystals will be reported