

## Enteric Pathogen Panel

**ALTERNATE NAMES**

Stool pathogens

**SPECIMEN TYPE**

Stool in Cary Blair transport media

**STORAGE**

Refrigerated

**ROUTINE / STAT T.A.T.**

24 hours

**REFERENCE RANGE**

Negative

**STABILITY**

Cairy Blair refrigerated: 48 hours

**SPECIAL INSTRUCTIONS**

Test includes:  
Campylobacter  
Salmonella  
Shigella  
Vibrio  
Yersinia  
Shigatoxin 1 and 2  
Norovirus  
Rotavirus