

MRSA SKIN + TISSUE

ALTERNATE NAMES

MRSA WOUND

SPECIMEN TYPE

SPECIAL RED SWAB
AVAILABLE FROM THE
LABORATORY

STORAGE

ROOM TEMP 24 HOURS
2-8 °C FOR 5 DAYS.

ROUTINE / STAT T.A.T.

1 DAY/3 HOURS

REFERENCE RANGE

NEGATIVE

STABILITY

ROOM TEMP 24 HOURS
2-8 °C FOR 5 DAYS

SPECIAL INSTRUCTIONS

THIS SWAB CAN ALSO BE
USED FOR AEROBIC CULTURE