

PSA

ALTERNATE NAMES

Prostate specific antigen

SPECIMEN TYPE

Serum

STORAGE

Refrigerate

ROUTINE / STAT T.A.T.

24 hours

REFERENCE RANGE

PSA Screen: ≤ 4.0 ng/ml
Post-cancer treatment: ≤ 0.1 ng/ml

STABILITY

3 hours room temp
24 hours refrigerated
5 months frozen

SPECIAL INSTRUCTIONS

Avoid hemolysis
Avoid turbidity