

RESPIRATORY PATHOGEN PANEL BY PCR

ALTERNATE NAMES

SPECIMEN TYPE

Nasal aspirates /washes or nasopharyngeal swabs collected and placed in universal transport medium.

STORAGE

2-8°C.

STABILITY

72 hours

ROUTINE / STAT T.A.T.

3 Hours

REFERENCE RANGE

Adenovirus: Negative
Human Metapneumovirus: Negative
Influenza A: Negative
Influenza A/H1: Negative
Influenza A/H3: Negative
Influenza B: Negative
Parainfluenza 1: Negative
Parainfluenza 2: Negative
Parainfluenza 3: Negative
Parainfluenza 4: Negative
RSV A: Negative
RSV B: Negative
Rhinovirus: Negative
B. parapertussis/bronchiseptica: Negative
B. homesii: Negative
B. pertussis: Negative

